From Diagnosis to Progress Of COVID-19



1.) Novel Coronavirus Detection Kit

Features

- 2019-nCoV's specific RNA detection
- Compatible with various RT-PCR device
- Manufactured in a KGMP facility
- Outstanding sensitivity: Excellent sensitivity based on high RNA/DNA affinity
- o High efficiency: Excellent PCR efficiency with specialized buffer system
- Quick Responsiveness: Improved Processivity and Elongation Rate

About Equipment Compatibility

- o BioRAD CFX96, CFX 96 Touch
- o ABI 7500 FAST
- Can be used with most RT-PCR Instruments







□ One-step Real-time RT-PCR Kit



KIT COMPONENTS	
Component Parts	Vol. (50rxn)
2* One Step RT-PCR Mix	650 ul (2ea)
Rnase Free Water	1ml
5* nCoV primer Mix-1 (RdRn gene)	250ul
5* nCoV primer Mix-2 (E gene)	250ul
Positive Control (P.C)	100ul





